

**NORTHWEST COMMUNITY EMS SYSTEM**  
**Non-Transport – Drug & Supply/Equipment List**  
**Last revised: 02-27-25 Rev.**

- EMS agencies shall assign appropriate personnel to inventory ambulances **daily** at shift change to ensure complete par levels, intact packaging, current dates, and good working order. All controlled substances must be viewed and counted daily per policy.
- The EMS MD or designees will do random unannounced non-transport inspections to measure compliance with these standards.
- All EMS products exchanged at hospitals must be **LATEX- FREE**. All non-exchange items must be latex-free unless a waiver has been granted and a latex-containing kit is maintained. Contain latex: Do NOT use without covering equipment or patient: BP cuffs, stethoscopes, Nellcor pulse oximeter.

Min.	ITEM
<b>IDPH Required BLS Supplies (see IDPH Non-Transport Inspection Form &amp; Section 515.825)</b>	
2	Adhesive tape rolls
1 ea.	Oropharyngeal airways: 5/100, 4/90, 3/80, 2/70, 1/60, 0/50
1 ea.	Nasopharyngeal airways: 12 Fr/3, 14 Fr/3.5, 16 Fr/4, 18 Fr/4.5, 20 Fr/5, 22 Fr/5.5, 24 Fr/6, 26 Fr/6.5, 28 Fr/7, 30 Fr/7.5
3	Water-soluble lubricant (individual packets)
2	Triangular bandages or slings
4	Roller bandages, self-adhering: 4" width
10	Sterile gauze pads: 4x4"
2	Vaseline gauze: 3x8"   Preferred alternative: Vented chest seal
2	Trauma/Universal dressings
1	Bandage shears
1	Emergency blanket (Mylar accepted)
1 ea.	Blood pressure cuff w/gauge: Adult, child, infant
1	Burn sheet
1 ea.	C-collars: Adult, child, infant
2	Cold packs
2	Hot packs
1	Flashlight/Pen light
1	OB kit w/ head cover and thermal absorbent blanket – aluminum foil or appropriate heat reflective material
1	Oxygen cylinder w/ min 350 liter capacity: min 1200psi with O2 key attached
1	Dial flowmeter/regulator for 15 liters per minute with delivery tubing
1 ea.	Oxygen masks: Adult, child, infant
1	Adult squeeze bag (BVM) with adult mask
1	Child squeeze bag (BVM) with child, infant and neonate mask
2	Splinting devices
1	Sterile saline or water solution (1000mL) in plastic bottle or bag
1	Stethoscope
1	Manually operated suction device
1	Defibrillator/AED w/ adult and pediatric pads
5	Run report forms
1	Glucometer (Microdot), test strips (bottle of 50), high/low test solution and at least 2 lancets
1	Res-Q-Pod
2	Tourniquet: Disposable Combat Application Tourniquet (CAT), TMT or approved equivalent
	Personal protection equipment (PPE) incl. nonporous disposable gloves, isolation bags, face/eye protection and gowns
	Communication equipment to reliably contact hospital
	ILS/ALS system approved equipment (medication storage box, airway equipment, monitor/defibrillator)

Min.	ITEM
<b>Medications/Supplies BLS: Above plus the following:</b>	
4	Acetaminophen PO chewable tablet 160mg or liquid equivalent
1	Albuterol 2.5mg/3mL (0.083%)
1	Ipratropium 0.5mg/2.5mL (0.02%)
1	Handheld nebulizer (HHN)
1	In-line nebulizer kit
4	ASA chewable 81mg
2	Diphenhydramine PO 50mg
2	Epinephrine 1:1 1mg/1mL
1	Glucagon 1mg
2	Naloxone 2mg/2mL
1	Nitroglycerin bottle 0.4mg
2	Ondansetron ODT 4mg
2	Medication atomization devices (MAD)
<b>Optional BLS</b>	
1	Automated CPR device: Defibtech Lifeline ARM, Stryker LUCAS 2/3, ROSC-U
<b>Medications/Supplies – ILS: Above plus the following:</b>	
2	Normal saline bags 1000mL
1 ea.	IV catheter (BD Nexiva): 14g, 16g, 18g, 20g
2	IV tubing: Macro drip 10, 15 or 20 gtts/mL
1	IV tubing: Micro drip 60 gtts/mL
2	IV protective dressings
4	CHG/IPA skin preps
1	Sharps container
1	Diphenhydramine IV 50mg/1mL
1	Dextrose 10% 25g/250mL
1	Ondansetron IV 4mg
<b>ILS: Selective drugs if vehicle has cardiac monitor with viewable ECG tracing, ETCO2 and SPO2</b>	
1	Adenosine 6mg/2mL
1	Amiodarone 150mg/3mL
1	Atropine sulfate 1mg/10mL
1	Midazolam 10mg/2mL
1	Fentanyl 100mcg/2mL
1	Lidocaine 2% 100mg/5mL
<b>Optional ILS</b>	
2	Disposable IV start pack with CHG/IPA skin prep
<b>Medications/Supplies – ALS: Above plus the following:</b>	
1	Adenosine 6mg/2mL
2	Amiodarone 150mg/3mL
2	Atropine sulfate 1mg/10mL
3	Epinephrine 1:10 1mg/10mL
1	Fentanyl 100mcg/2mL
1	Ketamine 500mg/10mL
1	Midazolam 10mg/2mL
1	Verapamil 5mg/2mL
1	10g catheter, 3" or 3.25" long for needle cric & pleural decompression
1 ea.	I-Gel: Size 5, Size 4, Size 3, Size 2.5, Size 2, Size 1.5

Min.	ITEM
1	Choking kit: (If King Vision/ProVu not carried: 1 standard laryngoscope handle and size 2 disposable blade); Magill or curved uterine forceps
1 ea.	Syringes: 1mL, 3mL, 10-12mL, 20mL with needles
1	Cardiac monitor/defibrillator: Equipped with adult and pediatric defibrillation pads capable of external pacing, SPO2, ETCO2 and 12-lead ECG
20	ECG electrodes
<b>Optional ALS</b>	
1 ea.	EZ-IO driver: 45mm (Yellow) needle, 25mm (Blue) needle, 15mm (Pink) needle with stabilizer kit
1	Sapphire infusion pump with current NWCEMSS Master Drug Library and 2 AP403 infusion sets
1	Flexicare ProVu Reusable 3.5" Video Monitor Display, charging cable and charging dock
2	ProVu Single Use Video Laryngoscope Handle with Hyper Blade Size 3.5
1	ProVu Single-Use Rigid Stylet Size Large
2 ea.	Endotracheal tubes (cuffed): 6 or 6.5mm, 7 or 7.5mm
1	ETT Holder
1	Battery operated suction unit with DuCanto catheters

**Mass Violence Bags – Issued per Region IX** (replacements responsibility of agency or System exchange if used)

#	Item	#	Item	#	Item
1	Maxpedition bag	2	Sharpie markers	4	Rescue blankets
1	Trauma shears	4	Lubricant (water soluble)	4	CAT or TMT tourniquet
4	28 Fr. NPA	4	Hyfin chest seal (more/agency expense)	3	Olaes bandage
2	Celox bandage (prefer <i>CeloxRapid Z fold</i> )	4	*Decompression Needles (strongly prefer minimum 10-12 g; 3" long needles) (Not issued for BLS bags)		

Recommend the availability of body armor (Helmet and vests) for those agencies participating in RTF teams.

## Vendor Contact Information:

<p><b>Bound tree: Sam Vollrath</b> Southern WI/Northern IL   Mobile: (608) 381-8172 <a href="mailto:Sam.Vollrath@Boundtree.com">Sam.Vollrath@Boundtree.com</a></p> <p><b>Angela Rashid</b> Account Manager, IL Phone: 614-760-5101   Fax: 866-561-1589 <a href="mailto:angela.rashid@boundtree.com">angela.rashid@boundtree.com</a></p>	<p><b>Eitan Group North America (Sapphire Pump)</b> <b>L.J. Relle, BBA, NRP, FP-C, CCP-C (Education)</b> <a href="mailto:lj@pacificbiomedical.com">lj@pacificbiomedical.com</a>   C:504-382-9638   <a href="http://www.eitanmedical.com">www.eitanmedical.com</a></p> <p><b>Purchases: Med Alliance: Steve Fina</b> 630-247-0770 <a href="mailto:stevefina@yahoo.com">stevefina@yahoo.com</a></p>
<p><b>Henry Schein EMS: Bill Beetschen</b> (847) 366-7989 (24/7)   <a href="mailto:bill.beetschen@henryschein.com">bill.beetschen@henryschein.com</a></p>	<p><b>Philips: Bill Carlson</b>   Cell: 224-230-9860 <a href="mailto:William.Carlson@Philips.Com">William.Carlson@Philips.Com</a></p>
<p><b>Tim Aebly: Account Executive DRE Medical</b> T: (877) 321-8194 / <a href="http://dremed.com">dremed.com</a> Carries certified refurbished Zoll &amp; PhysioControl monitors</p>	<p><b>Bill Cunningham: Penn Care, Inc.</b> Phone: (800) 392-7233 X 135 <a href="mailto:bill@penncare.net">bill@penncare.net</a> Best price on Bolin chest seals, LifePak 15 ECG Paper</p>
<p><b>Flexicare ProVu Laryngoscope</b> Caleb Waugh <a href="mailto:caleb.waugh@flexicare.com">caleb.waugh@flexicare.com</a> C 714-614-1234   T 949-450-9999</p>	<p><b>SafeBVM™ Adult Sotair device</b> Haris Shekhani <a href="mailto:haris@safebvm.com">haris@safebvm.com</a> 314-973-6606</p>
<p><b>Innovative Trauma Care: iTClamp</b>; 3463 Magic Drive, Suite 120 San Antonio, TX, 78229   1-855-774-4526 <a href="mailto:maggie@innovativetrauma.com">maggie@innovativetrauma.com</a> 727-483-4336 <a href="https://tactical-medicine.com/products/itclamp-hemorrhage-control-device">https://tactical-medicine.com/products/itclamp-hemorrhage-control-device</a></p>	
<p><b>Alix Webber - MDS - Vascular Access</b> E: <a href="mailto:alix.webber@bd.com">alix.webber@bd.com</a>   C: 847-239-2813</p>	<p><b>Safeguard medical: System pricing Celox and TMT tourniquet</b> <a href="mailto:jenniferro@safeguardmedical.com">jenniferro@safeguardmedical.com</a></p>
<p><b>Mercury Medical: Stormy Grinnell</b> Territory Mgr IL &amp; Iowa   Cell – (630) 432.1466 <a href="http://www.mercury-med.com">www.mercury-med.com</a></p>	<p><b>Teleflex (EZ-IO): Tyler Provost</b>   M: 216.894.7708 <a href="mailto:tyler.provost@teleflex.com">tyler.provost@teleflex.com</a>   Toll-free: 866-246-6990 North Chicago Sales Representative, Emergency Medicine</p>
<p><b>Quantum EMS Solutions</b> 516-321-9494   <a href="http://www.quantum-ems.com">www.quantum-ems.com</a>   <a href="mailto:sales@quantum-ems.com">sales@quantum-ems.com</a></p>	<p><b>David Dufek XDcuff®</b> P. (386) 275 -8361 <a href="mailto:david@xdcuff.com">david@xdcuff.com</a> Contact for special pricing</p>
<p><b>PhysioControl/Stryker Emerg Care: Steven Schmidt</b> C: (847) 636 0802 <a href="mailto:steven.schmidt@stryker.com">steven.schmidt@stryker.com</a></p> <p><b>Gianna Scarpetti</b> Sales Mgr – Chicagoland <b>Stryker: Patient Resuscitation and Treatment</b> C: 760 715 7025   <a href="http://www.lucas-cpr.com">www.lucas-cpr.com</a> <a href="mailto:gianna.scarpetti@stryker.com">gianna.scarpetti@stryker.com</a></p> <p><b>Savannah Johnson – Sales Rep</b> (616)-258-9426   <a href="mailto:savannah.johnson1@stryker.com">savannah.johnson1@stryker.com</a></p>	<p><b>Zoll</b> : Caroline Guibord Senior Account Executive EMS Cell: (773) 425-0710   <a href="mailto:cguibord@zoll.com">cguibord@zoll.com</a></p> <p><b>Danny Deutsch; Territory Manager – Ventilation</b> Phone: (248) 915-0122   <a href="mailto:Daniel.Deutsch@zoll.com">Daniel.Deutsch@zoll.com</a> Order supplies - <a href="https://www.zollwebstore.com/">https://www.zollwebstore.com/</a> Technical support – (800)348-9011</p>
<b>Disinfection Products</b>	
<p><b>AeroClave, LLC</b> Toll free: (800) 788-9119; Phone: (407) 788-3300 <a href="http://www.aeroclave.com">www.aeroclave.com</a></p>	<p><b>TOMI™ U.S. HEADQUARTERS (Steramist)</b> Phone: 800.525.1698   Fax: 301.710.5222 TOMIMIST.COM   <a href="mailto:info@tomimist.com">info@tomimist.com</a></p>
<p><b>Annette Uda: Aerapy, LLC</b> (UV light device   (630) 377-0348 x110; M: (630) 802-5244   <a href="http://www.Aerapy.com">www.Aerapy.com</a>   Meticulous cleaning and disinfecting of all surfaces per guidelines is required before using UV products.</p>	